



AIR FILTRATION BY LIFEGLOBAL®

Coda® Xtra Inline® Filters

The VOC removers – with a proven track record



Coda Xtra Inline filters

Unique VOC removing filter technology with a proven track record

- **AldaSorb®** removes over 98% of all aldehydes, formaldehydes and acetaldehydes from incoming gas lines
- **AldaSorb®** is safer and more effective in removing aldehydes than potassium permanganate
- Contain a HEPA filter to remove particulate matter
- Studies have proven that Coda® Inline® Filters are effective removing the majority of CACs (common air contaminants) and VOCs (volatile organic compounds)
- FDA cleared, CE certified and manufactured in a clean room facility in the USA
- Patented technology; ISO 13485 certified



Coda Xtra Inline Blue

Coconut shell-based activate carbon - change every 3 months (can be used up to 6 months); 20 times filtering material compared to standard



Coda Xtra Inline Green

Mixture of activated carbon, **AldaSorb®** - change every 3 months (can be used up to 6 months); 25 times filtering material compared to standard



Coda Regular Inline Purple

Improved activated carbon - change every 2 months (can be used up to 3 months); smaller footprint for labs with limited space; 10 times the filtering material compared to standard

| Product Code | Product Name | Quantity |
|--------------|-------------------------------|------------------|
| CXGR-001 | Coda Xtra Inline Green | One per gas line |
| CXBL-001 | Coda Xtra Inline Blue | One per gas line |
| CIPR-001 | Coda Regular Inline Purple | One per gas line |
| GILC-035 | Coda Inline Filter Connectors | One per gas line |
| GMAG-045 | Coda Inline Filter Bracket | One per gas line |