


4-Well GPS® Dish

BEST DECISION

less risk.

less worries.

less time.



embryo GPS®



4-Well GPS®



µDrop GPS®




Universal GPS®



38Special GPS®



Mini GPS®



embryo corral®

4-Well GPS® Dish

Key Features and Big Benefits

- 4 large wells with GPS bottoms for easy viewing and manipulation
- GPS bottom for embryo settling in the center of the wells for quick location of the embryo
- Wells are designed to reduce shadows and have concise focal positioning
- 4 large wells will hold up to 300 microliters of media and oil
- 4 smaller outer wells hold up to 100 microliter of media and oil
- Wide well design allows quick viewing and oil overlay of each well as well as the entire dish
- No pyrogenics, clean medical grade
- CE Certified and ISO 13485 certified
- 1-cell MEA and < 0.03 LAL Tested

Applications

- Holding oocytes prior to ICSI
- Large Volume Culture and Group Culture
- Assisted Hatching using laser
- Denuding before ICSI
- Cryopreservation and Thawing

IVFonline.com®

 **LifeGlobal®**
GROUP

T: 1-800-720-6375

F: 1-519-826-6947

Intl.: 001-519-826-5800

sales@LifeGlobal.com

www.LifeGlobalGroup.com

©2017 LifeGlobal® Group, LLC

Dish_4-Well-2017-08-16